Measuring potential differences

Recorder	N	/Janager			
Skeptic					
1 Setting up the voltmeter					
You will use a digital multimeter to me measure different things, you need to se	•		he meter can		
Plug the red probe lead into the rightmo	ost socket, labelle	d VΩHz.			v
Plug the black lead into the socket labell	led COM.				
Set the dial to Volts DC (the symbol with 200m. This means that the largest poten is 200 mV. Don't use scales other than V	tial difference tha	at can be measured	on this scale		
You may have expectations about the o your actual measurements, not your exp		ollowing measuren	nents. Record	10 A mA CO	M VΩHz
2 Measuring ∆ <i>V</i>				009	φ
Record the voltmeter readings you see w	hen you do the fo	ollowing:		Black	Red
a) Touch the voltmeter leads together:				wire	wire
b) Get a connecting wire (with alligator of to the red lead, and the other end to the				it. Clip one end	of the wire
c) Take a single battery out of the battery lead to the other end. What is the reading			end of the batter	ry marked +, an	d the black
d) The left justified "1" or "-1" is the volucial. Change the dial to point to the 2V	,			tmeter is set to	too small a
e) Reverse the leads of the voltmeter, so	the black lead is t	couched to the posi	itive end of the	battery:	
f) To get a positive reading, which voltme	eter lead should t	ouch the negative	end of the batte	ery?	
g) Use a ruler to measure the length of tery? Show your calculation:	the battery. What	is the magnitude of	of the average e	lectric field ins	ide the bat-
h) Take the other battery out of the hole tive end of the other. Measure the poten				nd of one touch	es the posi-
			6)-	+ 1)+	

i) What scale did you	have to use in part (h) to get a readin	ng other than "error"?	
j) Turn one of the bat	teries around and repeat the measure	ement:	+ + +
k) In the kit are two k	inds of light bulbs: short round light	bulbs, and long light	bulbs.
ONE connecting wire	ry, ONE ROUND light bulb, and from the kit. (Do not use a socket Make the light bulb light. Draw a circuit:		
CHECK WITH AND THEY DID THIS.	OTHER GROUP TO SEE HOW		
two connecting wires ROUND light bulb, c Measure and record to	back into the battery holder. Using s (with alligator clips), one battery, connect the circuit shown at right, so the following potential differences and the magnitude of the electric field instance.	o that the bulb is lit. In the length of each	A B C
-	Keep the leads in the same relative person as should be positive, others negative	-	
l) From A to B across	the battery		
ΔV =	L =	<u>a</u> =	
m) From B to C (use t	the most sensitive scale on the voltme	eter)	
ΔV =	L =	<u>2</u> =	
n) From C to D (the le	ength of the actual filament, uncoiled	d, is 1 cm)	
ΔV =	L = 0.01m (uncoiled)	$ \vec{E} =$	
Is the filament in equ	uilibrium? How do you know?		
o) From D to A (use the	he most sensitive scale on the voltmet	ter)	
ΔV =	L =	<u>`</u> =	
p) What is the round surements from parts	trip potential difference from A to B (l) through o):	to C to D back to A? S	Show your calculation, using your mea
g) Take the bulb and	socket out of the circuit. Measure the	e potential difference	across the bulb.
$\Delta V =$		librium?	
	Make sure that everyone in the Check with another group, the		